

Product datasheet for **RC201321L4V**

BOD1 (NM_138369) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	BOD1 (NM_138369) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BOD1
Synonyms:	FAM44B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138369
ORF Size:	555 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201321).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138369.1
RefSeq Size:	1934 bp
RefSeq ORF:	558 bp
Locus ID:	91272
UniProt ID:	Q96IK1
Cytogenetics:	5q35.2
MW:	19.2 kDa



[View online »](#)

Gene Summary:

Required for proper chromosome biorientation through the detection or correction of syntelic attachments in mitotic spindles.[UniProtKB/Swiss-Prot Function]