

## Product datasheet for **RC201330L1V**

### DCP1A (NM\_018403) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	DCP1A (NM_018403) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DCP1A
Synonyms:	HSA275986; Nbla00360; SMAD4IP1; SMIF
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_018403
ORF Size:	1746 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201330).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018403.3</a>
RefSeq Size:	6027 bp
RefSeq ORF:	1749 bp
Locus ID:	55802
UniProt ID:	<a href="#">Q9NPI6</a>
Cytogenetics:	3p21.1
Protein Families:	Transcription Factors
Protein Pathways:	RNA degradation



[View online »](#)

**MW:** 63.3 kDa

**Gene Summary:** Decapping is a key step in general and regulated mRNA decay. The protein encoded by this gene is a decapping enzyme. This protein and another decapping enzyme form a decapping complex, which interacts with the nonsense-mediated decay factor hUpf1 and may be recruited to mRNAs containing premature termination codons. This protein also participates in the TGF-beta signaling pathway. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Feb 2014]