

Product datasheet for **RC201338A1V**

Human PEX11C (PEX11G) (NM_080662) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human PEX11C (PEX11G) (NM_080662) AAV Particle
Tag: Myc-DDK
Symbol: PEX11C
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC201338 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGTCGCTGAGCGGCCTGGCGTGGCGCTGGAGTCGTACAGGGGCCGGGACCGCCTGATCCGAGTGC
TGGGGTACTGCTGCCAGCTGGTTGGTGGAGTTCTGGTTGAACAGTGTCCCGCCAGGTCCGAAGTGGGGAC
ACGTCTGTTGGTGGTGTCCACCAACTCAGCCACTGCAGGACCATTCTGCGACTCTTTGATGACCTGGCC
ATGTTTGTCTACACTAAGCAATATGGCCTGGGGGCACAGGAGGAGGACGCCTTTGTCCGCTGTGTCTCCG
TCCTAGGGAACCTGGCTGACCAGCTCTACTACCCTGTGAGCACGTGGCCTGGGCGGCTGATGCCCGGGT
CCTCCACGTGGACTCTTCTCGGTGGTGGACGCTGAGTACAACCTGTGGGCCCTCTCTCTGCTCCTGGGG
GTTGCCAGGTCCCTGTGGATGCTGTGAACTGAGACAGAGGCTGCGGAGCCCCACGGCGCCCTTCACCA
GCCCCGTGCCCGGGGCAAGCGGAGGGCCATGGAGGCGCAGATGCAGTCGGAGGCGCTGTCACTTCTCAG
CAACCTGGCCGACCTGGCCAACGCCGTGCACTGGCTGCCCCGGGGCGTGTGTGGCCGGCCGCTTCCCG
CCGTGGCTAGTGGGCTCATGGCACCATTCTCTCAATCCTCAGCATGTACCAGGCGGCCCGGGCCGGCG
GCCAGGCCGAGGCCACTACCCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201338 protein sequence
Red=Cloning site Green=Tags(s)

MASLSGLASALESYRGRDRLIRVLGYCCQLVGGVLVEQCPARSEVGTLLVVSTQLSHCRTILRLFDDLA
MFVYTKQYGLGAQEEDAFVRCVSVLGNLADQLYYPCEHVAWAADARVLHVDSSRWWTLSLTLWALSLLLG
VARSLWMLLKLRQRLRSPTAPFTSPLPRGKRRAMEAQMQSEALSLLSNLADLANAVHWLPRGVLWAGRFP
PWLVLGMGTISSILSMYQAARAGGQAEATTP

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_080662
ORF Size:	723 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_080662.2</u>
RefSeq Size:	1086 bp
RefSeq ORF:	726 bp
Locus ID:	92960
UniProt ID:	<u>Q96HA9</u>
Cytogenetics:	19p13.2
MW:	26.6 kDa