

Product datasheet for **RC201339L1V**

TNIP2 (NM_024309) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TNIP2 (NM_024309) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TNIP2
Synonyms:	ABIN2; FLIP1; KLIP
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_024309
ORF Size:	1287 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201339).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024309.2 , NP_077285.2
RefSeq Size:	1993 bp
RefSeq ORF:	1290 bp
Locus ID:	79155
UniProt ID:	Q8NFZ5
Cytogenetics:	4p16.3
MW:	48.7 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein which acts as an inhibitor of NFkappaB activation. The encoded protein is also involved in MAP/ERK signaling pathway in specific cell types. It may be involved in apoptosis of endothelial cells. Alternative splicing results in multiple transcript variants. A pseudogene related to this gene is located on the X chromosome.[provided by RefSeq, May 2014]