

Product datasheet for **RC201340A1V**

Human DOHH (NM_031304) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human DOHH (NM_031304) AAV Particle
Tag:	Myc-DDK
Symbol:	DOHH
Synonyms:	hDOHH; HLRC1
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201340 representing NM_031304 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGTGACGGAGCAGGAGGTGGATGCCATCGGGCAGACGCTGGTGGACCCCAAGCAGCCCTGCAGGCC
GCTTCCGGGCGCTGTTACGCTGCGTGGGCTCGGCGGCCAGGCGCCATTGCATGGATCAGCCAGGCCTT
CGATGACGATTCGCCCTGCTCAAGCAGAGCTGGCCTACTGCCTGGGCCAGATGCAGGATGCCCGGCC
ATCCCCATGCTGGTGGACGTGCTGAAGACACCCGTGAGGAGCCCATGGTGCCTCATGAGGCAGGGGAGG
CCCTGGGGGCCATCGGGACCCGGAAGTTCTGGAGATCCTGAAGCAGTATTCCTCGACCCCGTCATCGA
GGTGGCCGAGACCTGCCAGCTGGCCGTGCGCAGGCTGGAGTGGCTGCAGCAGCACGGCGGGGAGCCGGCG
GCGGGACCCTACCTCTCCGTGGACCCTGCCCGCCGGCTGAGGAGCGTGACGTGGGGCGCCTGCGGGAGG
CGCTGCTGGATGAGTCCCGGCCGCTCTTCGAGCGATACCGCGCCATGTTGCGCCTGCGCAACGCGGGAGG
CGAGGAGGCCGCCCTGGCGCTGGCCGAGGGTCTGCACTGTGGGAGCGCCCTTCCGCCACGAGGTCGGC
TACGTCTGGGACAGCTGCAGCACGAGGCGGCGGTGCCCGAGCTGGCGGCCCGCCTGGCCCGATGCACCG
AGAACCCCATGGTGGGCACGAGTGCGCGGAGGCCCTGGCGGCCATTGCCCGGCCCGCTGCCTGGCCGC
GCTGCAGGCTCACGCGGACGACCCAGAGCGCGTGGTGCCTGAGAGCTGCGAGGTGGCTCTGGACATGTAT
GAGCACGAGACCGGGCGGCCTTCCAGTACGCGGACGGCCTGGAGCAGCTGCGCGGGGCCCTCTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201340 representing NM_031304 Red=Cloning site Green=Tags(s) MVTEQEYDAIGQTLVDPKQPLQARFRALFTLRGLGGPGAIAWISQAFDDDSALLKHELAYCLGQMQDARA IPMLVDVLQDTRQEPMPVRHEAGEALGAIGDPEVLEILKQYSSDPVIEAETCQLAVRRLEWLQQHGGEPA AGPYLSYDPAPPAEERDVGRLEALLDESRPLFERYRAMFALRNAGGEEAALALAEGLHCGSALFRHEVG YVLGQLQHEAAVPQLAAALARCTENPMVRHECAEALGAIARPACLAALQAHADDPERVVRESCEVALDMY EHETGRAFYADGLEQLRGAPS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_031304
ORF Size:	906 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_031304.2</u>
RefSeq Size:	1772 bp
RefSeq ORF:	909 bp
Locus ID:	83475
UniProt ID:	<u>Q9BU89</u>
Cytogenetics:	19p13.3
MW:	32.7 kDa