

Product datasheet for **RC201362L4V**

HMBS (NM_000190) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	HMBS (NM_000190) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HMBS
Synonyms:	PBG-D; PBGD; PORC; UPS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000190
ORF Size:	1083 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201362).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000190.3
RefSeq Size:	1526 bp
RefSeq ORF:	1086 bp
Locus ID:	3145
UniProt ID:	P08397
Cytogenetics:	11q23.3
Domains:	Porphobil_deam
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Metabolic pathways, Porphyrin and chlorophyll metabolism

MW: 39.3 kDa

Gene Summary: This gene encodes a member of the hydroxymethylbilane synthase superfamily. The encoded protein is the third enzyme of the heme biosynthetic pathway and catalyzes the head to tail condensation of four porphobilinogen molecules into the linear hydroxymethylbilane. Mutations in this gene are associated with the autosomal dominant disease acute intermittent porphyria. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008]