

Product datasheet for **RC201371L4V**

UBAP2L (NM_014847) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	UBAP2L (NM_014847) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UBAP2L
Synonyms:	NICE-4; NICE4
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_014847
ORF Size:	3261 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201371).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_014847.2 , NP_055662.2
RefSeq Size:	4006 bp
RefSeq ORF:	3264 bp
Locus ID:	9898
UniProt ID:	Q14157
Cytogenetics:	1q21.3
Domains:	UBA
Protein Families:	Druggable Genome



[View online »](#)

MW: 114.5 kDa

Gene Summary: Plays an important role in the activity of long-term repopulating hematopoietic stem cells (LT-HSCs).[UniProtKB/Swiss-Prot Function]