

Product datasheet for RC201379A1V

Human SNURF (NM_005678) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SNURF (NM_005678) AAV Particle
Tag: Myc-DDK
Symbol: SNURF
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC201379 representing NM_005678
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGAGCGGGCAAGGGATCGCTTACACCTGAGACGAACTACAGAACAGCACGTACCAGAGGTGGAAGTCC
 AAGTCAAACGCAGAAGGACTGCCTCACTGAGCAACCAAGAGTGTCAAGTTGTACCCGAGGCGTTCTCAGCA
 GCAGCAAGTACCTGTGGTGGATTCCAGGCTGAAGTGAAGCAGGCATTCTAGCTGAGACACCAAGAGGT
 GGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201379 representing NM_005678
 Red=Cloning site Green=Tags(s)
 MERARDRLHLRRTTEQHVPEVEVQVKRRRTASLSNQECQLYPRRSQQQVPVDFQAE LRQAFLAETPRG
 G

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_005678
ORF Size: 213 bp
Buffer: PBS with 0.001% Pluronic F68


[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_005678.3
RefSeq Size:	1326 bp
RefSeq ORF:	216 bp
Locus ID:	8926
UniProt ID:	Q9Y675
Cytogenetics:	15q11.2
MW:	8.2 kDa