

Product datasheet for **RC201382L3V**

DPH2 (NM_001384) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DPH2 (NM_001384) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DPH2
Synonyms:	DPH2L2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001384
ORF Size:	1467 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201382).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001384.4
RefSeq Size:	2513 bp
RefSeq ORF:	1470 bp
Locus ID:	1802
UniProt ID:	Q9BQC3
Cytogenetics:	1p34.1
Domains:	Diphthamide_syn
MW:	52.1 kDa



[View online »](#)

Gene Summary:

This gene is one of two human genes similar to the yeast gene dph2. The yeast gene was identified by its ability to complement a diphthamide mutant strain, and thus probably functions in diphthamide biosynthesis. Diphthamide is a post-translationally modified histidine residue present in elongation factor 2 (EF2) that is the target of diphtheria toxin ADP-ribosylation. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2016]