

Product datasheet for **RC201388A1V**

Human ELAC1 (NM_018696) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human ELAC1 (NM_018696) AAV Particle
Tag:	Myc-DDK
Symbol:	ELAC1
Synonyms:	D29
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201388 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCTATGGATGTGACATTCTGGGGACGGGTGCAGCATACCCATCTCCAACCCGGGGTGCCTCTGCTG
TGGTCCTTCGGTGTGAAGGCGAGTGTGGCTCTTTGACTGTGGGGAGGGAACACAGACACAGCTTATGAA
AAGCCAACCTAAAGCAGGGAGAATTACCAAGATCTTCATCACACACCTTCATGGAGACCATTTCTTTGGC
CTTCTGGGCTCCTCTGCACAATCAGCCTGCAGAGTGGCTCCATGGTGTCCAAACAGCCTATTGAAATCT
ATGGCCCTGTAGGGCTTCGGGACTTTATCTGGCGAACCATGGAATCTCTCACACGGAGCTGGTCTTCCA
TTATGTGGTTCATGAACTGGTTCCTACAGCAGATCAATGTCCTGCAGAAGAACTAAAAGAATTTGCGCAT
GTGAATAGAGCAGACAGTCTCCCAAAGAGGAACAAGGAAGAACTATCCTGTTAGACTCAGAAGAAAAC
CATACCTTCTGTTTGATGATGAACAATTTGTTGTAAGCAATTCGCCTCTTTCACAGAATTCCTCATT
TGGGTTTTCTAGTCGTGAAAAGAAACGCCAGGTAACTCAATGCACAGAACTTAAAGACCTTGGTGT
CCACCAGGTCTGCCTATGGGAAGCTGAAAAATGGAATTTCTGTTGTTCTGAAAATGGGTTACAATTT
CTCCCAAGATGTCTTAAAAAGCCTATTGTTGGAAGAAAAATCTGCATATTGGGTGACTGCTCTGGGT
TGTGGGTGATGGAGGAGTAAACTGTGCTTTGAAGCAGACCTGTTGATCCAGCAAGCAACCTGGATGAT
GCCAGATGGACAAAGCAAAGGAGCATGGCCACAGCACACCACAGATGGCAGCAACATTTGCAAAGTTGT
GCCGTGCAAGAGGCTGGTTCGACTCACTTCAGTCAGAGGTACAAACCAGTTGCCTTGGCCAGAGAAGG
AGAAACAGATGGCATTGCAGAACTAAAAAGCAAGCTGAATCAGTGTTAGATCTCCAAGAAGTGACTCTA
GCAGAAGATTTTATGGTGATAAGCATTCCAATCAAGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201388 protein sequence Red=Cloning site Green=Tags(s) MSMDVTFGLGTGAAYPSPTRGASAVVLRCEGECWLFDCGEGTQTQLMKSQLKAGRITKIFITHLHGDHFFG LPGLLCTISLQSGSMVSKQPIEYGPVGLRDFIWRTMELSHTLVFHYVVHEL VPTADQCPAEELKEFAH VNRADSPPKKEEQGRTILLDSEENS YLLFDDEQFVVKAFRLFHRIPSFGFSVVEKKRPGKLN AQLKDLGV PPGPAYGKLKNGISVVL ENGVTISPQDVLKKPIVGRKICILGDCSGVVG DGGVKLCFEADLLIHEATLDD AQMDKAKEHGHSTPQAATFAKLCRAKRLVLTHFSQRYKPVALAREGETDGI AELKKQAESVLDLQEVTL AEDFMVISIPIKK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_018696
ORF Size:	1089 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_018696.2
RefSeq Size:	2233 bp
RefSeq ORF:	1092 bp
Locus ID:	55520
UniProt ID:	Q9H777
Cytogenetics:	18q21.2
MW:	40 kDa