

Product datasheet for **RC201389A1V**

Human YIPF2 (NM_024029) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human YIPF2 (NM_024029) AAV Particle
Tag:	Myc-DDK
Symbol:	YIPF2
Synonyms:	FinGER2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201389 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCATCGGCCGACGAGCTGACCTTCCATGAATTCGAGGAGGCCACTAATCTTCTGGCTGACACCCAG
ATGCAGCCACCACCAGCAGAAGCGATCAGCTGACCCACAAGGGCACGTGGCTGTGGCCGTGGGCTCAGG
TGGCAGCTATGGAGCCGAGGATGAGGTGGAGGAGGAGTGAAGGCCGCGCTCCTGCAGGAGCAGCAG
CAGCAGCAGCAGCCGGGATTCTGGACCTTCAGCTACTATCAGAGCTTCTTTGACGTGGACACCTCACAGG
TCCTGGACCGGATCAAAGGCTCACTGCTGCCCCGGCCTGGCCACAACCTTTGTGCGGCACCATCTGCGGAA
TCGGCCGGATCTGTATGGCCCCCTTCTGGATCTGTGCCACGTTGGCCTTTGTCTGCGCCGTCAGTGCAAC
CTGACGCTGGTGTGGCCAGAGGAGGGACCCCTCCATCCACTACAGCCCCAGTTCACAAGGTGACCG
TGGCAGGCATCAGCATCTACTGCTATGCGTGGCTGGTGGCCCTGGCCCTGTGGGGCTTCTGCGGTGGCG
CAAGGGTGTCCAGGAGCGCATGGGGCCCTACACCTTCTGGAGACTGTGTGCATCTACGGCTACTCCCTC
TTTGTCTTCATCCCCATGGTGGTCTGTGGCTCATCCCTGTGCCTTGGCTGCAGTGGCTCTTTGGGGCGC
TGGCCCTGGGCTGTGAGCCGCCGGGCTGGTATTCACCTCTGGCCCGTGGTCCGTGAGGACACCAAGCT
GGTGGCCACAGTGTGCTGTCCGTGGTGTGCTGCTCCAGCCCTCCTGGCCATGGGCTGTAAGTTGTAC
TTCTTCAGTCGCTGCCTCCGAGAACGTGGCTCCTCCTCCCAATCACATCTCTGCCCTCAAACATCG
CGCTGTCCCTACCTTGCCGAGTCCCTGGCCCCCTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201389 protein sequence Red=Cloning site Green=Tags(s) MASADELTFHEFEEATNLLADTPDAATTSRSDQLTPQGHVAVAVGSGGSYGAEDEVEEESDKAALLQEQQ QQQQPGFWTFSYYQSFFDVDTSQVLDRIKGSLLPRPGHNFVRHHLRNRPDLYGPFWICATLAFVLAVTGN LTLVLAQRDPSTHYSPQFHKVTVAGISYCYAWLVPLALWGFLRWKGVQERMGPYTFLETVCYGYSL FVFIPMVVLWLIPVPWLQWLFALALGLSAAGLVFTLWPVVREDTRLVATVLLSVVLLHALLAMGCKLY FFQSLPPENVAPPPQITSLPSNIALSPTLPQSLAPS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_024029
ORF Size:	948 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_024029.3 , NP_076934.1
RefSeq Size:	1554 bp
RefSeq ORF:	951 bp
Locus ID:	78992
UniProt ID:	Q9BWQ6
Cytogenetics:	19p13.2
MW:	35.2 kDa