

## Product datasheet for RC201393A1V

### Human ELOB (NM\_007108) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human ELOB (NM\_007108) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** ELOB  
**Synonyms:** SIII; TCEB2  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC201393 ORF sequence  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGATCGCC**

ATGGACGTGTTCTCATGATCCGGCGCCACAAGACCACCATCTTCACGGACGCCAAGGAGTCCAGCACGG  
 TGTTCGAACTGAAGCGCATCGTCGAGGGCATCCTCAAGCGGCCTCTGACGAGCAGCGGCTGTACAAGGA  
 TGACCAACTCTTGATGATGGCAAGACTGGGCGAGTGTGGCTTCACCAGTCAAACAGCACGGCCACAG  
 GCCCCAGCCACAGTGGGGCTGGCCTTCGGGCAGATGACACCTTTGAGGCCCTGTGCATCGAGCCGTTTT  
 CCAGCCCGCCAGAGCTGCCGATGTGATGAAGCCCCAGGACTCGGGAAGCAGTGCCAATGAACAAGCCGT  
 GCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC201393 protein sequence  
 Red=Cloning site Green=Tags(s)

MDVFLMIRRHKTTIFTDAKESSTVFELKRIVEGILKRPPDEQRLYKDDQLLDDGKTLGECGFTSQTPARQ  
 APATVGLAFRADDTFEALCIEPFSSPELPDVMKPDQSGSSANEQAVQ

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_007108



[View online »](#)

<b>ORF Size:</b>	354 bp
<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_007108.2</a></u>
<b>RefSeq Size:</b>	1009 bp
<b>RefSeq ORF:</b>	357 bp
<b>Locus ID:</b>	6923
<b>UniProt ID:</b>	<u><a href="#">Q15370</a></u>
<b>Cytogenetics:</b>	16p13.3
<b>MW:</b>	13.1 kDa