

Product datasheet for **RC201400L4V**

FOXSI (NM_004118) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	FOXSI (NM_004118) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FOXSI
Synonyms:	FKHL18; FREAC10
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_004118
ORF Size:	990 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201400).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004118.3
RefSeq Size:	1353 bp
RefSeq ORF:	993 bp
Locus ID:	2307
UniProt ID:	O43638
Cytogenetics:	20q11.21
Protein Families:	Transcription Factors
MW:	35.4 kDa


[View online »](#)

Gene Summary:

The forkhead family of transcription factors belongs to the winged helix class of DNA-binding proteins. The protein encoded by this intronless gene contains a forkhead domain and is found predominantly in aorta and kidney. The function of the encoded protein is unknown. [provided by RefSeq, Jul 2008]