

Product datasheet for **RC201403L4V**

ABCF3 (NM_018358) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ABCF3 (NM_018358) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ABCF3
Synonyms:	EST201864
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018358
ORF Size:	2127 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201403).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018358.1
RefSeq Size:	2626 bp
RefSeq ORF:	2130 bp
Locus ID:	55324
UniProt ID:	Q9NUQ8
Cytogenetics:	3q27.1
Domains:	ABC_tran, AAA
Protein Families:	Druggable Genome



[View online »](#)

MW: 79.7 kDa

Gene Summary: This gene encodes a member of the ATP-binding cassette (ABC) transporter superfamily. ATP-binding cassette proteins transport various molecules across extra- and intracellular membranes. The protein encoded by this gene displays antiviral effect against flaviviruses. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2017]