

Product datasheet for RC201407

NAP1L4 (NM_005969) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NAP1L4 (NM_005969) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NAP1L4
Synonyms:	hNAP2; NAP1L4b; NAP2; NAP2L
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC201407 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGCAGATCACAGTTTTTCAGATGGGGTTCCTTCAGATTCGGTGAAGCTGCTAAAAATGCAAGTAACA
 CAGAAAAGCTCACAGATCAGGTGATGCAGAATCCTCGAGTTCTGGCAGCTTTACAGGAGCGACTTGACAA
 TGTCCCTCACACCCCTTCCAGCTACATCGAACTTTACCTAAAGCAGTAAAAAGAAGAATTAATGCATTG
 AAACAACCTTCAGGTGAGATGTGCTCACATAGAAGCCAAGTTCTATGAAGAGGTACATGACTTGGAAGAA
 AGTATGCAGCGCTATACCAGCCTCTCTTTGACAAGAGAAGAGAATTTATCACCGCGCATGTTGAACCAAC
 AGATGCGGAATCGGAATGGCACAGTGAAATGAAGAGGAAGAGAAATTGGCTGGAGACATGAAAAGTAAA
 GTAGTCGTACAGAAAAAGCAGCGGCAACGGCTGAAGAGCCAGATCCCAAAGGAATTCAGAGTTCTGGT
 TTACCATCTTCAGAAATGTGGACATGCTGAGTGAATTAGTCCAGGAATATGATGAACCAATCTTGAAACA
 CCTGCAGGATATTAAGTGAAATTTCTGACCCTGGACAGCCTATGTCTTTGTGTTAGAGTTCCACTTT
 GAACCAACGACTACTTTACCAACTCAGTCCTGACAAAAACCTACAAGATGAAATCAGAACCAGATAAGG
 CTGATCCCTTTTCCTTTGAAGTCCTGAGATTGTGGACTGTGACGGGTGTACTATTGACTGGAAGAAAGG
 AAAGAATGTTACTGTCAAAACCATCAAGAAAAAGCAGAAGCATAAGGGTCGAGGCACTGTTAGAACAAAT
 ACGAAACAAGTACCAATGAGTCCTTTTCAACTTCTTCAATCCATTGAAAGCATCCGGGGATGGAGAAT
 CACTGGATGAAGATTCTGAATTCACATTAGCCTCTGATTTTGAAATTGGACACTTTTCCGTGAGCGGAT
 AGTCCCGCGGGTGTGCTGTACTTCACTGGGGAGGCCATAGAAGATGATGACAATTTTGAAGAAGGTGAA
 GAAGGAGAAGAGGAGGAATTAGAAGGTGACGAGGAGGAGAAGACGAGGATGATGCGGAAATTAACCCCA
 AGGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RC201407 protein sequence
 Red=Cloning site Green=Tags(s)

MADHSFSDGVPSDSVEAAKNASNTEKLTDQVMQNPRVLAALQERLDNVPHTPSSYIETLPKAVKRRINAL
 KQLQVRCAHIEAKFYEEVHDLERKYAALYQPLFDKRREFITGDVEPTDAESEWHSNEEEEEKLAGDMKSK
 VVVTEKAAATAEPPDPKGIPEFWFTIFRNVDMLSLVQEYDEPILKHLQDIKVKFSDPGQPMSEVLEFHF
 EPNDYFTNSVLTKTYKMKSEPKADPFSFEGPEIVDCDGCTIDWKKGKNVTVKTIKKKQKHKGRGTVRTI
 TKQVPNESFFNFFNPLKASGDGSLDESEFTLASDFEIGHFFRERIVPRAVLYFTGEAIEDDDNFEEGE
 EEEEELEGEDEGEDDAEINPKV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6576_h03.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_005969

ORF Size: 1125 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.


Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_005969.4
RefSeq Size:	2564 bp
RefSeq ORF:	1128 bp
Locus ID:	4676
UniProt ID:	Q99733
Cytogenetics:	11p15.4
Domains:	NAP
MW:	42.8 kDa
Gene Summary:	This gene encodes a member of the nucleosome assembly protein (NAP) family which can interact with both core and linker histones. It can shuttle between the cytoplasm and nucleus, suggesting a role as a histone chaperone. This gene is one of several located near the imprinted gene domain of 11p15.5, an important tumor-suppressor gene region. Alterations in this region have been associated with the Beckwith-Wiedemann syndrome, Wilms tumor, rhabdomyosarcoma, adrenocortical carcinoma, and lung, ovarian, and breast cancer. [provided by RefSeq, Jul 2008]


Product images:



Circular map for RC201407



Western blot validation of overexpression lysate (Cat# [LY416953]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC201407 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified NAP1L4 protein (Cat# [TP301407]). The protein was produced from HEK293T cells transfected with NAP1L4 cDNA clone (Cat# RC201407) using MegaTran 2.0 (Cat# [TT210002]).