

Product datasheet for **RC201416L4V**

C17orf75 (NM_022344) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C17orf75 (NM_022344) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C17orf75
Synonyms:	NJMU-R1; SRI2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022344
ORF Size:	1188 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201416).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022344.2
RefSeq Size:	4599 bp
RefSeq ORF:	1191 bp
Locus ID:	64149
UniProt ID:	Q9HAS0
Cytogenetics:	17q11.2
Protein Families:	Druggable Genome
MW:	44.6 kDa



[View online »](#)

Gene Summary:

As component of the WDR11 complex acts together with TBC1D23 to facilitate the golgin-mediated capture of vesicles generated using AP-1 (PubMed:29426865). May have a role in spermatogenesis.[UniProtKB/Swiss-Prot Function]