

Product datasheet for **RC201449**

zinc finger protein 775 (ZNF775) (NM_173680) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	zinc finger protein 775 (ZNF775) (NM_173680) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	zinc finger protein 775
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC201449 ORF sequence, **codon optimized**.
 Due to the complexity of NM_173680, the ORF clone is codon optimized for mammalian Expression.
 The nucleotide sequence differs from the reference sequence, yet the amino acid sequence remains identical.

Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCC**GCGATCGCC**

ATGAGAGCGGACTTGCAGGGAATGGCACAGGAGCTGGGCTTGTTCATGAAGGTCAAACAAGAAAAGCCCGAACGCCTCTTGCAGACCTGGCGCCCCAGGCCATGCTGGTTGAGAAAGACAAAGAAAACATTTTTTCAGCAACAGAGGGTTGCCACCAAGACAGACGATGGGGAGGCCCGGGCCCTCGGTGGACAGGAGGAATCTGGAAGTCTAGATGGGCCCTCCAACAGAGCAAGACGCTGGTCTGGCCGGACGAGCTCCGGGCTCCGCTTCAGGCCCCCTGAGTCCTAGCCTTAGCAGCGGTGAAGGACACTTCGCTGCCTGGATTGCGGAAAACGCTTCTCTTGGTGGTCATCACTGAAAATACACCAACGCACACTCACACTGGAGAAAAGCCATATCTTTGTGGGAAGTGGGAAGTCTTTCTCTCAGAAGCCTAACCTGGCCCGACATCAGCGGCATCACACAGGGGAGCGGCCGTTTTGCTGTCCAGAAATGCGCACGGCGGTTTTCCCAAAAACAGCACCTGCTGAAGCACCAGAAGACACAGTCGGCAGCGCACTCCTGCCCTGAATGCGAACGGTGTCTTCGACACCAGGTCGGTCTTAGAATCCATCAGCGGGCCACGCAAGGGATAGGCAGGGATCAAGAGCTGGGCTGCACGAGCTCATCCAGGATGCTGCCGCACGCCGAGCTTCCCGCTCCAGCTGGCCCCCAGAGGGCGCCGGAGTGGGCTGGCTGGGACTGTGTCAAGGCTGGTGGGACAACCAGGCGCACGCCCGCTGTCCGGACCTGAGGGGCCAGGGAAACCGCCATTTATTTGCAACGAGTGTGGAAAGTCTTTACCTGGTGGTCCAGCCTGAACATTACCAACGCATTACATCCGGGAAAGCCATATGCTTGCCAGAATGTGGCCGAGGTTTTCCAGAAACAACTTGACCCGGCACCTCCGGAATCACACCGGAGAGAGGCCCCACCCCTGCCCTCACTGCGGCCGCGGTTTTCGGCAGAAGCAGCATCTCCTTAAGCACCTCCGCACCCATCTTCCAGGAGCCAGGCTGCGCCCTGTCCCTCATGTGGCAAGTCTGCGGCTCACGCGCAGATTGCGCGCCACCAGCGGGCACACGCGGTGGCTGAGCCGGCAGTTCGGCAGGAGCCAGGGGACCAGCCTCAGGCGGAAGCCATCCCCGGCTGGCTGCTAGGCCCGGTCTCTCAGAGATCTCCAGGCGCAAGAGATACCCTTTGGGACGCGGGCAGGCTGGCCTGGCCGGTCCCGGGAGCCTAGACAGTTTATATGTAACGAGTGTGGCAAATCCTTCTCCTGGTGGTCTGCTGACTATCCATCAGAGGATCCATACAGGAGAAAGCCATACCCTTGCCCGGAATGCGGAAGACGCTTCTCCAGAAAGCCAAATTTGACACGCCATCGCAGAAATCATACAGGAGAGCGACCATACTTGTGTCCCGCTTGTGGGAGGGGTTTTCCAGAAACAACATCTCTTGAAACACCAGAGGGTTCATCGGGCAGCACCCGCTTGCAGCCAAAAGAGGAAGCGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>RC201449 representing NM_173680
 Red=Cloning site Green=Tags(s)

MESGLAGNGTGAGLVMKVQEKPERLLQTLAPQAMLVEKDKENIFQQHRGLPPRQTMGRPRALGGQEESGSPRWAPTEQDAGLAGRAPGSASGPLSPSLSSGEGHFVCLDCGKRFSSWSSSLKIHQRTHTGKPYLPCGCKGKSFQKPNLARHQRHHTGERPFCCPECARRFQKQHLKHKQKTHSRPATHSCPECERCFRHQVGLRIHQRAHARDRQGSRAGLHEL IQDAAARRACRLQGPGRPEWAWLGLCQGWGQPGARAASVSGPEGPGEPRQFICNECGKSFTWSSLNIIHQRIHTGERPYACPECGRRFSQKPNLTRHLRNHTGERPHPCPHCGRFRQKHLLKHLRTHLPGAQAAPCPCGKSCRRAALRAHQRAHAVAEPVAPAGEPGDQPQAEIPLAARPRSSQRSPGARDTLWGRGQAGLAGPGEPRQFICNECGKSFSSWWSALTIHQRIHTGERPYPCPECGRRFSQKPNLTRHRNHTGERPYLCPACGRGFSQKQHLKHKQVHRAAPACSPKEEAR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_173680

ORF Size: 1611 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_173680.3](#), [NP_775951.2](#)

RefSeq Size: 2261 bp


RefSeq ORF: 1614 bp

Locus ID: 285971

UniProt ID: [Q96BV0](#)

Cytogenetics: 7q36.1
 MW: 59.8 kDa
 Gene Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC201449