

## Product datasheet for **RC201458A1V**

### Human HMGA1 (NM\_002131) AAV Particle

#### Product data:

**Product Type:** AAV Particles  
**Product Name:** Human HMGA1 (NM\_002131) AAV Particle  
**Tag:** Myc-DDK  
**Symbol:** HMGA1  
**Synonyms:** HMG-R; HMGA1A; HMG1Y  
**Mammalian Cell Selection:** None  
**Vector:** pAAV-AC-Myc-DDK (PS100089)  
**ORF Nucleotide Sequence:** >RC201458 ORF sequence  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAGTGAGTCGAGCTCGAAGTCCAGCCAGCCCTTGGCCTCCAAGCAGGAAAAGGACGGCACTGAGAAGC  
GGGGCCGGGGCAGGCCGCGCAAGCAGCCTCCGAAGGAGCCAGCGAAGTGCCAACACCTAAGAGACCTCG  
GGGCCGACCAAGGGAAGCAAAAACAAGGGTGTGCAAGACCCGGAAAACCACCACAACCTCCAGGAAGG  
AAACCAAGGGGCAGACCCAAAAAAGTGGAGAAGGAGGAAGAGGAGGGCATCTCGCAGGAGTCTCGGAGG  
AGGAGCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RC201458 protein sequence  
Red=Cloning site Green=Tags(s)

MSESSKSSQPLASKQEKDGTGTEKRGRGRPRKQPPKEPSEVPTPKRPRGRPKGSKNKGAATRKTTPGR  
KPRGRPKKLEKEEEEGISQESSEEEQ

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

**Species:** Human  
**Serotype:** AAV-2  
**ACCN:** NM\_002131  
**ORF Size:** 288 bp



[View online »](#)

<b>Buffer:</b>	PBS with 0.001% Pluronic F68
<b>Stability:</b>	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
<b>RefSeq:</b>	<u><a href="#">NM_002131.2</a></u>
<b>RefSeq Size:</b>	1998 bp
<b>RefSeq ORF:</b>	291 bp
<b>Locus ID:</b>	3159
<b>UniProt ID:</b>	<u><a href="#">P17096</a></u>
<b>Cytogenetics:</b>	6p21.31
<b>MW:</b>	10.7 kDa