

## Product datasheet for **RC201485L4V**

### **BXDC2 (BRIX1) (NM\_018321) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	BXDC2 (BRIX1) (NM_018321) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BXDC2
Synonyms:	BRIX; BXDC2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018321
ORF Size:	1059 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201485).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018321.3</a>
RefSeq Size:	1303 bp
RefSeq ORF:	1062 bp
Locus ID:	55299
UniProt ID:	<a href="#">Q8TDN6</a>
Cytogenetics:	5p13.2
Domains:	Brix
Protein Families:	Embryonic stem cells, Stem cell - Pluripotency



[View online »](#)

**MW:** 41.4 kDa

**Gene Summary:** Required for biogenesis of the 60S ribosomal subunit.[UniProtKB/Swiss-Prot Function]