

Product datasheet for RC201488L1

TMEM33 (NM_018126) Human Tagged Lenti ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	TMEM33 (NM_018126) Human Tagged Lenti ORF Clone
Tag:	Myc-DDK
Symbol:	TMEM33
Synonyms:	1600019D15Rik; Pom33; SHINC-3; SHINC3
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201488).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF.


ACCN:	NM_018126
ORF Size:	741 bp




[View online >](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_018126.1
RefSeq Size:	2684 bp
RefSeq ORF:	744 bp
Locus ID:	55161
UniProt ID:	P57088
Cytogenetics:	4p13
Domains:	UPF0121
Protein Families:	Transmembrane
MW:	27.8 kDa
Gene Summary:	Acts as a regulator of the tubular endoplasmic reticulum (ER) network. Suppresses the RTN3/4-induced formation of the ER tubules (PubMed:25612671). Positively regulates PERK-mediated and IRE1-mediated unfolded protein response signaling (PubMed:26268696). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for RC201488L1



Double digestion of RC201488L1 using SgfI and MluI