

## Product datasheet for **RC201492L3V**

### **RBM22 (NM\_018047) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RBM22 (NM_018047) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RBM22
Synonyms:	Cwc2; fSAP47; ZC3H16
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_018047
ORF Size:	1260 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201492).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_018047.1</a>
RefSeq Size:	2352 bp
RefSeq ORF:	1263 bp
Locus ID:	55696
UniProt ID:	<a href="#">Q9NW64</a>
Cytogenetics:	5q33.1
Domains:	RRM, zf-CCCH
Protein Pathways:	Spliceosome



[View online »](#)

**MW:** 46.9 kDa

**Gene Summary:** This gene encodes an RNA binding protein. The encoded protein may play a role in cell division and may be involved in pre-mRNA splicing. Related pseudogenes exist on chromosomes 6, 7, 9, 13, 16, 18, and X. [provided by RefSeq, Mar 2009]