

Product datasheet for **RC201503A1V**

Human TXNL4B (NM_017853) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human TXNL4B (NM_017853) AAV Particle
Tag: Myc-DDK
Symbol: TXNL4B
Synonyms: Dim2; DLP
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC201503 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGCTTCCTACTGCCCAAGCTGACTAGCAAAAAGGAAGTAGACCAGGCGATAAAAAGTACTGCTGAGA
AGGTGTTGGTTCTCAGGTTTGGGAGAGATGAAGATCCTGTCTGTCTGCAGCTAGATGATATTCTTTCTAA
GACCTCTTCTGACTTAAGTAAAAAGGCTGCTATATACCTGGTAGATGTGGACCAAAGTGCAGTTTATACA
CAGTATTTTGACATCAGTTATATTCCATCTACTGTCTTTTTCTCAATGGCAGCATATGAAAGTGGATT
ATGGATCTCCAGATCACACTAAGTTTGTGGGAAGCTTCAAACCAACAAGACTTCATAGATTTGATTGA
AGTAATCTATCGAGGAGCAATGAGGGGAAGCTTATTGTCCAAAGTCTATTGATCCCAAGAATATTCCC
AAATATGACCTTCTCTATCAAGACATT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201503 protein sequence
Red=Cloning site Green=Tags(s)

MSFLLPKLTSKKEVDQAIKSTAEKVLVLRFRDEDPVCLQLDDILSKTSSDLKMAAIYLVVDVDTAVYT
QYFDISYIPSTVFFNGQHMKVDYGPDPHTKFGVSKTKQDFIDLIEVIYRGAMRGKLIQSPIDPKNIP
KYDLLYQDI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online >](#)

ACCN:	NM_017853
ORF Size:	447 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_017853.1</u>
RefSeq Size:	2563 bp
RefSeq ORF:	450 bp
Locus ID:	54957
UniProt ID:	<u>Q9NX01</u>
Cytogenetics:	16q22.2
MW:	17 kDa