

## Product datasheet for **RC201521L1V**

### GLMN (NM\_053274) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	GLMN (NM_053274) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GLMN
Synonyms:	FAP; FAP48; FAP68; FKBPAP; GLML; GVM; VMGLOM
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_053274
ORF Size:	1782 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201521).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053274.2</a>
RefSeq Size:	2042 bp
RefSeq ORF:	1785 bp
Locus ID:	11146
UniProt ID:	<a href="#">Q92990</a>
Cytogenetics:	1p22.1
Protein Families:	Druggable Genome
MW:	68.2 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a phosphorylated protein that is a member of a Skp1-Cullin-F-box-like complex. The protein is essential for normal development of the vasculature and mutations in this gene have been associated with glomuvenous malformations, also called glomangiomas. Multiple splice variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2016]