

### Product datasheet for RC201522L1

#### 9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

## PGM3 (NM\_015599) Human Tagged Lenti ORF Clone

**Product data:** 

**Product Type:** Expression Plasmids

Product Name: PGM3 (NM 015599) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: PGM3

Synonyms: AGM1; IMD23; PAGM; PGM 3

Mammalian Cell None

Selection:

Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC201522).

Sequence:

**Restriction Sites:** Sgfl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF.

**ACCN:** NM\_015599

ORF Size: 1626 bp





#### PGM3 (NM\_015599) Human Tagged Lenti ORF Clone - RC201522L1

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** <u>NM 015599.1</u>

RefSeq Size:6107 bpRefSeq ORF:1629 bpLocus ID:5238

 UniProt ID:
 O95394

 Cytogenetics:
 6q14.1

Domains: PGM PMM, PGM PMM I

**Protein Pathways:** Amino sugar and nucleotide sugar metabolism

**MW:** 59.9 kDa


**Gene Summary:** This gene encodes a member of the phosphohexose mutase family. The encoded protein

mediates both glycogen formation and utilization by catalyzing the interconversion of glucose-1-phosphate and glucose-6-phosphate. A non-synonymous single nucleotide polymorphism in this gene may play a role in resistance to diabetic nephropathy and neuropathy. Alternatively spliced transcript variants encoding multiple isoforms have been


observed for this gene. [provided by RefSeq, Dec 2010]



# **Product images:**



Circular map for RC201522L1



Double digestion of RC201522L1 using Sgfl-Mlul