

Product datasheet for **RC201537L1V**

NFIL3 (NM_005384) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NFIL3 (NM_005384) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NFIL3
Synonyms:	E4BP4; IL3BP1; NF-IL3A; NFIL3A
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_005384
ORF Size:	1386 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201537).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005384.2
RefSeq Size:	2104 bp
RefSeq ORF:	1389 bp
Locus ID:	4783
UniProt ID:	Q16649
Cytogenetics:	9q22.31
Domains:	BRLZ
Protein Families:	Transcription Factors



[View online »](#)

MW: 51.5 kDa

Gene Summary: The protein encoded by this gene is a transcriptional regulator that binds as a homodimer to activating transcription factor (ATF) sites in many cellular and viral promoters. The encoded protein represses PER1 and PER2 expression and therefore plays a role in the regulation of circadian rhythm. Three transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Feb 2014]