

Product datasheet for **RC201543A1V**

Human N myc interactor (NMI) (NM_004688) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human N myc interactor (NMI) (NM_004688) AAV Particle
Tag:	Myc-DDK
Symbol:	N myc interactor
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201543 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGAAGCTGATAAAGATGACACACAACAAATCTTAAGGAGCATTGCCAGATGAATTTATAAAAGATG
AACAAATAAGGGACTAATTGATGAAATTACAAAGAAAAATTTCAACTAAAGAAGGAGATCCAAAAGCT
TGAAACGGAGTTACAAGAGGCTACCAAGAATTCCAGATTAAGAGGATATTCTGAAACAAAGATGAAA
TTCTTATCAGTTGAACTCCTGAGAATGACAGCCAGTTGTCAAATATCTCCTGTTCAATTTCAAGTGAGCT
CGAAAGTTCCTTATGAGATACAAAAGGACAAGCACTTATCACCTTTGAAAAAGAAGAAGTTGCTCAAAA
TGTGGTAAGCATGAGTAAACATCATGTACAGATAAAAGATGTAATCTGGAGGTTACGGCCAAGCCAGTT
CCATTAATTCAGGAGTCAGATTCAGGTTTATGTAGAAGTTTCTAAAATGAAAAATCAATGTTACTGAAA
TTCCTGACACATTGCGTGAAGATCAAATGAGAGACAACTAGAGCTGAGCTTTTCAAAGTCCCGAAATGG
AGGCGGAGAGGTGGACCGCTGGACTATGACAGACAGTCCGGGAGTGCAGTCATCACGTTTGTGGAGATT
GGAGTGGCTGACAAGATTTTGAAAAAGAAAGAAATACCCTCTTTATATAAATCAAACCTGCCATAGAGTTA
CTGTTTCTCCATACACAGAAATACACTTGAAAAAGTATCAGATATTTTCAGGAACATCTAAGAGGACAGT
GCTTCTGACAGGAATGGAAGGCATTCAAATGGATGAAGAAATTTGTGGAGGATTTAATTAACATTCACTTT
CAACGGGCAAAGAATGGAGGTGGAGAAGTAGATGTGGTCAAGTGTCTCTAGGTCAACCTCACATAGCAT
ACTTTGAAGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201543 protein sequence Red=Cloning site Green=Tags(s) MEADKDDTQQILKEHSPDEFIKDEQNKGLIDEITKKNIQLKKEIQKLETQLQEATKEFQIKEDIPETKMK FLSVETPENDSQLSNISCSFQVSSKVPYEQKGQALITFEKEEVAQNVSMSKHHVQIKDYNLEVTAKPV PLNSGVRFQVYVEVSKMKINVTEIPDTLREDQMRDKLELSFSKSRNGGGEVDRVDYDRQSGSAVITFVEI GVADKILKKKEYPLYINQTCRVTVPYTEIHLKKYQIFSGTSKRTVLLTGMEGIQMDEEIVEDLINIHF QRAKNGGGEVDVVKCSLGQPHIAYFEE TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_004688
ORF Size:	921 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_004688.1
RefSeq Size:	1501 bp
RefSeq ORF:	924 bp
Locus ID:	9111
UniProt ID:	Q13287
Cytogenetics:	2q23.3
MW:	35.1 kDa