

Product datasheet for **RC201563L2V**

EIF2S2 (NM_003908) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | EIF2S2 (NM_003908) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | EIF2S2 |
| Synonyms: | eIF-2-beta; EIF2; EIF2B; EIF2beta; PPP1R67 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_003908 |
| ORF Size: | 999 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC201563). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_003908.2 |
| RefSeq Size: | 2592 bp |
| RefSeq ORF: | 1002 bp |
| Locus ID: | 8894 |
| UniProt ID: | P20042 |
| Cytogenetics: | 20q11.22 |
| Domains: | eIF2B_5 |
| MW: | 38.4 kDa |


[View online »](#)

Gene Summary:

Eukaryotic translation initiation factor 2 (EIF-2) functions in the early steps of protein synthesis by forming a ternary complex with GTP and initiator tRNA and binding to a 40S ribosomal subunit. EIF-2 is composed of three subunits, alpha, beta, and gamma, with the protein encoded by this gene representing the beta subunit. The beta subunit catalyzes the exchange of GDP for GTP, which recycles the EIF-2 complex for another round of initiation. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2015]