

Product datasheet for RC201565A1V

Human CSRP2 (NM_001321) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CSRP2 (NM_001321) AAV Particle
Tag: Myc-DDK
Symbol: CSRP2
Synonyms: CRP2; LMO5; SmLIM
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC201565 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC

ATGCTGTCTGGGGAGGTGGAACAAGTGTGGGCCTGTGGGAGGACCGTGTACCACGCAGAAGAGGTGC
 AGTGTGATGGCAGGAGCTTCCACCGCTGCTGCTTTCTCTGCATGGTTGCAGGAAAAATTTAGATAGCAC
 AACAGTGGCAATTCACGATGAAGAGATCTACTGCAATCCTGCTACGGAAGAAGTATGGGCCAAAAGGC
 TACGGTTATGGCCAGGGCGCTGGCAGCCTTAACATGGACCGTGGCGAGAGGCTGGGCATCAAACCAGAGA
 GTGTTACGCCTCACAGGCCTACAACAAATCCAAACACTTCTAAATTTGCTCAGAAATATGGAGGTGCTGA
 GAAGTGTTCAGATGTGGGATTCTGTATATGCTGCCGAGAAGATAATTGGAGCTGGAAGCCCTGGCAC
 AAAAACTGTTTCCGATGTGCAAAGTGTGGGAAGAGTCTTGAATCAACAACCTCTGACTGAAAAAGAAGGTG
 AAATCTATTGTAAAGGATGCTATGCAAAGAACTTTGGGCCCAAGGGATTGGCTATGGCCAAGGAGCAGG
 GGCTCTTGTTTCATGCCAG

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201565 protein sequence
 Red=Cloning site Green=Tags(s)

MPVWGGGNKCGACGRTVYHAEVQCDGRSFHRCFLCMVCRKNLDSTTVAIHDEEYCKSCYKKGYPKG
 YGYGQGAGTLNMDRGERLGIKPESVQPHRPTTNPNTSKFAQKYGGAEKCSRCGDSVYAAEKIIGAGKPWH
 KNCFRCAKCGKSLESTLTKEGEYCKGKYAKNFGPKGFGYQGAGALVHAQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human


[View online »](#)

Serotype:	AAV-2
ACCN:	NM_001321
ORF Size:	579 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001321.1</u>
RefSeq Size:	940 bp
RefSeq ORF:	582 bp
Locus ID:	1466
UniProt ID:	<u>Q16527</u>
Cytogenetics:	12q21.2
MW:	21 kDa