

Product datasheet for **RC201569L2V**

Annexin A1 (ANXA1) (NM_000700) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | Annexin A1 (ANXA1) (NM_000700) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | Annexin A1 |
| Synonyms: | ANX1; LPC1 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_000700 |
| ORF Size: | 1038 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC201569). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_000700.1 |
| RefSeq Size: | 1399 bp |
| RefSeq ORF: | 1041 bp |
| Locus ID: | 301 |
| UniProt ID: | P04083 |
| Cytogenetics: | 9q21.13 |
| Domains: | annexin |
| Protein Families: | Druggable Genome |



[View online »](#)

MW: 38.7 kDa

Gene Summary: This gene encodes a membrane-localized protein that binds phospholipids. This protein inhibits phospholipase A2 and has anti-inflammatory activity. Loss of function or expression of this gene has been detected in multiple tumors. [provided by RefSeq, Dec 2014]