

Product datasheet for **RC201586L3V**

SOX15 (NM_006942) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SOX15 (NM_006942) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SOX15
Synonyms:	SOX20; SOX26; SOX27
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006942
ORF Size:	699 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201586).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006942.1
RefSeq Size:	1394 bp
RefSeq ORF:	702 bp
Locus ID:	6665
UniProt ID:	O60248
Cytogenetics:	17p13.1
Protein Families:	Transcription Factors
MW:	25.1 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the SOX (SRY-related HMG-box) family of transcription factors involved in the regulation of embryonic development and in the determination of the cell fate. The encoded protein may act as a transcriptional regulator after forming a protein complex with other proteins. [provided by RefSeq, Jul 2008]