

Product datasheet for **RC201596A1V**

Human SNAP23 (NM_003825) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human SNAP23 (NM_003825) AAV Particle
Tag: Myc-DDK
Symbol: SNAP23
Synonyms: HsT17016; SNAP-23; SNAP23A; SNAP23B
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC201596 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGATAATCTGTCATCAGAAGAAATCAACAGAGAGCTCACCAGATTACTGATGAGTCTCTGGAAAGTA
CGAGGAGAATCCTGGGTTTAGCCATTGAGTCTCAGGATGCAGGAATCAAGACCATCACTATGCTGGATGA
ACAAAAGGAACAACTAAACCGCATAGAAGAAGGCTTGGACCAATAAATAAGGACATGAGAGAGACAGAG
AAGACTTTAACAGAACTCAACAAATGCTGTGGCCTTTGTGTCTGCCCATGTAATAGAACAAGAACTTTG
AGTCTGGCAAGGCTTATAAGACAACATGGGGAGATGGTGGAGAAAACTCACCTTGCAATGTAGTATCTAA
ACAGCCAGGCCCGGTGACAAATGGTCAGCTTCAGCAACCAACAACGGGAGCAGCCAGTGGTGGATACATT
AAACGCATAACTAATGATGCCAGAGAAGATGAAATGGAAGAGAACCTGACTCAAGTGGGCAGTATCCTGG
GAAATCTAAAAGACATGGCCCTGAACATAGGCAATGAGATTGATGCTCAAAATCCACAAATAAACGAAT
CACAGACAAGGCTGACACCAACAGAGATCGTATTGATATTGCCAATGCCAGAGCAAGAAACTCATTGAC
AGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201596 protein sequence
Red=Cloning site Green=Tags(s)

MDNLSSEIIQRAHQITDESLESTRILGLAIESQDAGIKTITMLDEQKEQLNRIEEGLDQINKDMRETE
KTLTELNKCCGLCVCPCNRTKNFESGKAYKTTWGDGGENSPCNVVSQPGPVTNGQLQQPTTGAASGGYI
KRITNDAREDEMEENLTQVGSILGNLKDMLNIGNEIDAQNPQIKRITDKADTNDRIDIANARAKKLID
S

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_003825
ORF Size:	633 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_003825.2
RefSeq Size:	2650 bp
RefSeq ORF:	636 bp
Locus ID:	8773
UniProt ID:	O00161
Cytogenetics:	15q15.1-q15.2
MW:	23.4 kDa