

Product datasheet for **RC201597L3V**

MAPK6 (NM_002748) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	MAPK6 (NM_002748) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MAPK6
Synonyms:	ERK3; HsT17250; p97MAPK; PRKM6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002748
ORF Size:	2163 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201597).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002748.2
RefSeq Size:	4225 bp
RefSeq ORF:	2166 bp
Locus ID:	5597
UniProt ID:	Q16659
Cytogenetics:	15q21.2
Domains:	pkinese, TyrKc, S_TKc
Protein Families:	Druggable Genome, ES Cell Differentiation/IPS, Protein Kinase



[View online »](#)

MW: 82.7 kDa

Gene Summary: The protein encoded by this gene is a member of the Ser/Thr protein kinase family, and is most closely related to mitogen-activated protein kinases (MAP kinases). MAP kinases also known as extracellular signal-regulated kinases (ERKs), are activated through protein phosphorylation cascades and act as integration points for multiple biochemical signals. This kinase is localized in the nucleus, and has been reported to be activated in fibroblasts upon treatment with serum or phorbol esters. [provided by RefSeq, Jul 2008]