

## Product datasheet for RC201605L1

#### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

### Proteasome 20S beta 6 (PSMB6) (NM\_002798) Human Tagged Lenti ORF Clone

**Product data:** 

**Product Type: Expression Plasmids** 

**Product Name:** Proteasome 20S beta 6 (PSMB6) (NM\_002798) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: Proteasome 20S beta 6

DELTA; LMPY; Y Synonyms:

**Mammalian Cell** 

Selection:

Vector: pLenti-C-Myc-DDK (PS100064) E. coli Selection:

None

Chloramphenicol (34 ug/mL) The ORF insert of this clone is exactly the same as(RC201605).

Sequence:

**ORF Nucleotide** 

Sgfl-Mlul

**Restriction Sites: Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF.

ACCN: NM\_002798

**ORF Size:** 717 bp



#### Proteasome 20S beta 6 (PSMB6) (NM\_002798) Human Tagged Lenti ORF Clone - RC201605L1

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** <u>NM 002798.1</u>

 RefSeq Size:
 869 bp

 RefSeq ORF:
 720 bp

 Locus ID:
 5694

 UniProt ID:
 P28072

 Cytogenetics:
 17p13.2

**Domains:** proteasome

**Protein Families:** Druggable Genome, Protease

Protein Pathways: Proteasome MW: 25.4 kDa


**Gene Summary:** The proteasome is a multicatalytic proteinase complex with a highly ordered ring-shaped 20S

core structure. The core structure is composed of 4 rings of 28 non-identical subunits; 2 rings are composed of 7 alpha subunits and 2 rings are composed of 7 beta subunits. Proteasomes are distributed throughout eukaryotic cells at a high concentration and cleave peptides in an ATP/ubiquitin-dependent process in a non-lysosomal pathway. The encoded protein is a member of the proteasome B-type family, also known as the T1B family, and is a 20S core beta subunit in the proteasome. Alternatively spliced transcript variants encoding multiple


isoforms have been observed for this gene. [provided by RefSeq, Jul 2012]



# **Product images:**



Circular map for RC201605L1



Double digestion of RC201605L1 using Sgfl and Mlul