

Product datasheet for **RC201614L3V**

GBA (NM_001005741) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GBA (NM_001005741) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GBA
Synonyms:	GBA1; GCB; GLUC
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001005741
ORF Size:	1608 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201614).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001005741.1
RefSeq Size:	2583 bp
RefSeq ORF:	1611 bp
Locus ID:	2629
UniProt ID:	P04062
Cytogenetics:	1q22
Protein Families:	Druggable Genome
Protein Pathways:	Lysosome, Metabolic pathways, Other glycan degradation, Sphingolipid metabolism



[View online »](#)

MW: 59.7 kDa

Gene Summary: This gene encodes a lysosomal membrane protein that cleaves the beta-glucosidic linkage of glycosylceramide, an intermediate in glycolipid metabolism. Mutations in this gene cause Gaucher disease, a lysosomal storage disease characterized by an accumulation of glucocerebrosides. A related pseudogene is approximately 12 kb downstream of this gene on chromosome 1. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2010]