

Product datasheet for RC201629L3V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Beclin 1 (BECN1) (NM_003766) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: Beclin 1 (BECN1) (NM_003766) Human Tagged ORF Clone Lentiviral Particle

Symbol: Beclin 1

Synonyms: ATG6; beclin1; VPS30

Mammalian Cell

Selection:

Puromycin

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

Tag: Myc-DDK
ACCN: NM 003766

ORF Size: 1350 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC201629).

Sequence:
OTI Disclaimer:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 003766.2

 RefSeq Size:
 2144 bp

 RefSeq ORF:
 1353 bp

 Locus ID:
 8678

 UniProt ID:
 Q14457

Cytogenetics: 17q21.31

Domains: APG6

Protein Families: Druggable Genome





Beclin 1 (BECN1) (NM_003766) Human Tagged ORF Clone Lentiviral Particle - RC201629L3V

Protein Pathways: Regulation of autophagy

MW: 51.7 kDa

Gene Summary: This gene encodes a protein that regulates autophagy, a catabolic process of degradation

induced by starvation. The encoded protein is a component of the phosphatidylinositol-3-kinase (PI3K) complex which mediates vesicle-trafficking processes. This protein is thought to play a role in multiple cellular processes, including tumorigenesis, neurodegeneration and apoptosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep

2015]