

Product datasheet for RC201631

RPL15 (NM_002948) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: RPL15 (NM_002948) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: RPL15
Synonyms: DBA12; EC45; L15; RPL10; RPLY10; RPYL10
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC201631 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGGTGCATACAAGTACATCCAGGAGCTATGGAGAAAGAAGCAGTCTGATGTCATGCGCTTTCTTCTGA
 GGGTCCGCTGCTGGCAGTACCGCCAGCTCTCTGCTCTCCACAGGGCTCCCCGCCCCACCCGGCCTGATAA
 AGCGCGCCGACTGGGCTACAAGGCCAAGCAAGGTTACGTTATATATAGGATTCGTGTTCCCGTGGTGCC
 CGAAAACGCCAGTTCCTAAGGGTGCAACTACGGCAAGCCTGTCCATCATGGTGTAAACCAGCTAAAGT
 TTGCTCGAAGCCTTCAGTCCGTTGCAGAGGAGCGAGCTGGACGCCACTGTGGGGCTCTGAGAGTCCTGAA
 TTCTTACTGGGTTGGTGAAGATTCCACATACAAATTTTTTGGAGTTATCCTCATTGATCCATTCCATAAA
 GCTATCAGAAGAAATCCTGACACCCAGTGGATCACCAAACAGTCCACAAGCACAGGGAGATGCGTGGGC
 TGACATCTGCAGGCCGAAAGAGCCGTGGCCTTGGAAAGGGCCACAAGTCCCACCACACTATTGGTGGCTC
 TCGCCGGGCGACTTGAGAAAGGCCAATACTCTCCAGTCCACCGTTACCGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201631 protein sequence
 Red=Cloning site Green=Tags(s)

MGAYKYIQELWRKKQSDVMRFLLRVRCWQYRQLSALHRAPRPTRPDKARRLGYKAKQGYVIYRIRVRRGG
 RKRVPVPGATYGKPVHHGVNQLKFARSLQSVAEERAGRHCALRVLNSYWVGEDSTYKFFEVILIDPFHK
 AIRRNPDTQWITKPVHKHREMRGLTSAGRKSRGLGKGKSHHTIGGSRAAWRRRNTLQLHRYR


TRTRPLEQKLISEEDLAANDILDYKDDDDKV



Chromatograms: https://cdn.origene.com/chromatograms/mk6323_e02.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_002948

ORF Size: 612 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_002948.5](#)

RefSeq Size: 2350 bp

RefSeq ORF: 615 bp

Locus ID: 6138

UniProt ID: [P61313](#)

Cytogenetics: 3p24.2


Domains: Ribosomal_L15e

Protein Pathways: Ribosome


MW: 24.1 kDa

Gene Summary: Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of four RNA species and approximately 80 structurally distinct proteins. This gene encodes a member of the L15E family of ribosomal proteins and a component of the 60S subunit. This gene shares sequence similarity with the yeast ribosomal protein YL10 gene. Elevated expression of this gene has been observed in esophageal tumors and gastric cancer tissues, and deletion of this gene has been observed in a Diamond-Blackfan anemia (DBA) patient. As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome. [provided by RefSeq, Mar 2017]

Product images:



Circular map for RC201631



Western blot validation of overexpression lysate (Cat# [LY418995]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC201631 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).