

Product datasheet for RC201645L1V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

PARK7 (NM_007262) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: PARK7 (NM_007262) Human Tagged ORF Clone Lentiviral Particle

Symbol: PARK7

Synonyms: DJ-1; DJ1; GATD2; HEL-S-67p

Mammalian Cell

Selection:

None

Vector: pLenti-C-Myc-DDK (PS100064)

Tag: Myc-DDK
ACCN: NM 007262

ORF Size: 567 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC201645).

Sequence:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This

naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 007262.3, NP 009193.2

 RefSeq Size:
 979 bp

 RefSeq ORF:
 570 bp

 Locus ID:
 11315

 UniProt ID:
 Q99497

 Cytogenetics:
 1p36.23

Domains: DJ-1_PfpI

Protein Families: Druggable Genome, Protease





PARK7 (NM_007262) Human Tagged ORF Clone Lentiviral Particle - RC201645L1V

Protein Pathways: Parkinson's disease

MW: 19.9 kDa

Gene Summary: The product of this gene belongs to the peptidase C56 family of proteins. It acts as a positive

regulator of androgen receptor-dependent transcription. It may also function as a redox-sensitive chaperone, as a sensor for oxidative stress, and it apparently protects neurons against oxidative stress and cell death. Defects in this gene are the cause of autosomal recessive early-onset Parkinson disease 7. Two transcript variants encoding the same protein

have been identified for this gene. [provided by RefSeq, Jul 2008]