

Product datasheet for **RC201655A1V**

Human RFC3 (NM_002915) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human RFC3 (NM_002915) AAV Particle
Tag:	Myc-DDK
Symbol:	RFC3
Synonyms:	RFC38
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201655 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGAGCCTCTGGGTGGACAAGTATCGGCCCTGCTCCTTGGGACGGCTGGACTATCACAAGGAGCAGGCGG
CCCAGCTGCGGAACCTGGTGCAGTGTGGTGACTTTCCTCATCTGTTAGTGTACGGACCATCAGGTGCTGG
AAAAAAGACAAGAATTATGTGTATTCTACGTGAACCTTTATGGTGTGGAGTGGAAAAATTGAGAATTGAA
CATCAGACCATCACAACCTCCATCTAAAAAATTGAAATTAGCACCATTGCAAGTAACCTACCACCTTG
AAGTTAATCCTAGTGATGCTGGAAATAGTGACCGAGTAGTCATTGAGGAGATGTTGAAAACAGTGGCACA
ATCACAACAACCTTGAACAACCTCTCAAGGGATTTTAAAGTGGTATTATTGACAGAAAGTTGACAAACTC
ACCAAAGATGCTCAGCATGCCTTGCGAAGAACCATGGAAAAATATATGTCTACCTGCAGATTGATCTTGT
GCTGCAATTCTACATCTAAAGTGATCCACCTATTCGTAGTAGGTGCTTGGCGGTTCTGTGCTGCTCC
CAGCATTGAAGATATTTGCCACGTGTTATCTACTGTGTGAAGAAGGAAGGTCTGAATCTTCCTTCACAA
CTGGCTCATAGACTTGACAGAGAAGTCTGTAGAAATCTCAGAAAAGCCCTGCTTATGTGTGAAGCCTGCA
GAGTGCAACAATATCCTTTACTGCAGATCAAGAAATCCCTGAGACAGATTGGGAGGTGTATCTGAGGGA
GACTGCAAAATGCTATTGTGAGTCAAGAACTCCACAAAGGCTCCTTGAAGTTCGTGGAAGGCTGTATGAG
CTTCTAACTATTGATTCCTCCTGAGATAATAATGAAGGGCTTCTCTCAGAACTGTTACATAATTGTG
ATGGACAACTGAAAGGGGAGGTGGCACAATGGCAGCTTACTATGAGCATCGTCTACAGCTGGGTAGCAA
AGCCATTTATCACTTGAAGCGTTTGTGGCCAAATTCATGGCACTTTATAAGAAGTTCATGGAGGATGGA
TTGGAAGGCATGATGTTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201655 protein sequence Red=Cloning site Green=Tags(s) MSLVWDKYRPCSLGRLDYHKEQAAQLRNLVQCGDFPHLLVYGPSGAGKKTRIMCILREL YGVGVEKLRIE HQTITTPSKKKIEISTIASNYHLEVNPSDAGNSDRVVIQEMLKTVAQSQQLETNSQRDFKVVLLTEVDKL TKDAQHALRRTMEKYMSTCRLILCCNSTSKVIPIRSRCLAVRVPAPSIEDICHVLSTVCKKEGLNLP SQ LAHRLAEKSCRNLRKALLMCEACRVQQYPFTADQEIPETDWEVYLRETANAIVSQQTPQRLLEVRGRLYE LLTHCIPPEIIMKGLLSELLHNCDGQLKGEVAQMAAYYEHRLQLGSKAIYHLEAFVAKFMALYKKFMEDG LEGMMF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_002915
ORF Size:	1068 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_002915.3
RefSeq Size:	2396 bp
RefSeq ORF:	1071 bp
Locus ID:	5983
UniProt ID:	P40938
Cytogenetics:	13q13.2
MW:	40.6 kDa