

Product datasheet for **RC201663L4V**

CEP72 (NM_018140) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CEP72 (NM_018140) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CEP72
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018140
ORF Size:	1941 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201663).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018140.3
RefSeq Size:	2427 bp
RefSeq ORF:	1944 bp
Locus ID:	55722
UniProt ID:	Q9P209
Cytogenetics:	5p15.33
MW:	71.7 kDa
Gene Summary:	The product of this gene is a member of the leucine-rich-repeat (LRR) superfamily of proteins. The protein is localized to the centrosome, a non-membraneous organelle that functions as the major microtubule-organizing center in animal cells. [provided by RefSeq, Jul 2008]



[View online »](#)