

Product datasheet for **RC201704A1V**

Human NUCKS1 (NM_022731) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human NUCKS1 (NM_022731) AAV Particle
Tag:	Myc-DDK
Symbol:	NUCKS1
Synonyms:	JC7; NUCKS
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201704 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCGCGGCCTGTCAGAAATAGGAAGTTGTTGATTACTCACAGTTTCAGGAATCTGATGATGCAGATG
AAGATTATGGAAGAGATTCGGGCCCTCCCACTAAGAAAATTCGATCATCTCCCGAGAAGCTAAAAATAA
GAGGCGATCTGGAAAGAATTCACAGGAAGATAGTGAGGACTCAGAAGACAAAGATGTGAAGACCAAGAAG
GATGATTCTCACTCAGCAGAGGATAGTGAAGATGAAAAAGAAGATCATAAAATGTGCGCAACAACGGC
AGGCGGCATCTAAGCAGCTTCTAAACAGAGAGAGATGCTCATGGAAGATGTGGGCAGTGAGGAAGAACA
AGAAGAGGAGGATGAGGCACCATTCAGGAGAAAGATTCCGGCAGCGATGAAGATTTCTAATGGAAGAT
GATGACGATAGTGACTATGGCAGTTGAAAAAGAAAAAAGATGGTTAAGAAGTCCAAACCTGAAA
GAAAAGAAAAGAAAATGCCAAACCCAGACTAAAGGCTACAGTGACGCCAAGTCCAGTGAAAGGCAAAGG
GAAAGTGGGTGCGCCACAGCTTCAAAGGCATCAAAGGAAAAGACTCCTTCTCCCAAAGAAGAAGATGAG
GAACCGGAAAGCCCGCCAGAAAAGAAAACATCTACAAGCCCCCACCAGAAATCTGGGGATGAAGGGT
CTGAAGATGAAGCCCCTTCTGGGGAGGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201704 protein sequence Red=Cloning site Green=Tags(s) MSRPVRNRKVVDSYQFQESDDADEDYGRDSGPPTKKIRSSPREAKNRRSGKNSQEDSEDSKDVKTKK DDSHSAEDSEDEKEDHKNVRQQRQAASKAASKQREMLMEDVGSEEEQEEDEAPFQEKDSGSDFLMED DDSDYGSKKKNKMKVKKSKPERKEKKMPKRLKATVTPSPVKGKGKVGRRPTASKASKEKTPSPKEEDE EPESPPEKKTSTSPPEKSGDEGSEDEAPSGED TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_022731
ORF Size:	729 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_022731.2</u>
RefSeq Size:	6471 bp
RefSeq ORF:	732 bp
Locus ID:	64710
UniProt ID:	<u>Q9H1E3</u>
Cytogenetics:	1q32.1
MW:	27.3 kDa