

Product datasheet for **RC201713L1V**

SNRP70 (SNRNP70) (NM_003089) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SNRP70 (SNRNP70) (NM_003089) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SNRP70 |
| Synonyms: | RNPU1Z; RPU1; Snp1; SNRP70; U1-70K; U1AP; U1RNP; U170K |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-Myc-DDK (PS100064) |
| Tag: | Myc-DDK |
| ACCN: | NM_003089 |
| ORF Size: | 1311 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC201713). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_003089.4 |
| RefSeq Size: | 2008 bp |
| RefSeq ORF: | 1314 bp |
| Locus ID: | 6625 |
| UniProt ID: | P08621 |
| Cytogenetics: | 19q13.33 |
| Protein Families: | Druggable Genome |
| Protein Pathways: | Spliceosome |



[View online »](#)

MW: 51.6 kDa

Gene Summary: Component of the spliceosomal U1 snRNP, which is essential for recognition of the pre-mRNA 5' splice-site and the subsequent assembly of the spliceosome (PubMed:19325628, PubMed:25555158). SNRNP70 binds to the loop I region of U1-snRNA (PubMed:2467746, PubMed:19325628, PubMed:25555158).[UniProtKB/Swiss-Prot Function]