

Product datasheet for **RC201772A1V**

Human HEAB (CLP1) (NM_006831) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human HEAB (CLP1) (NM_006831) AAV Particle
Tag:	Myc-DDK
Symbol:	HEAB
Synonyms:	hClp1; HEAB
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201772 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGGAGAAGAGGCTAATGATGACAAGAAGCCAACCACTAAATTTGAACTAGAGCGAGAAACAGAACTTC
GCTTTGAGGTGGAGGCATCTCAGTCAGTTGAGTTGTTGACTGGCATGGCAGAGATCTTTGGCAC
AGAGCTGACCCGAAACAAGAAATTCACCTTTGATGCTGGTGCCAAAGGTGGCTGTTTTCACTTGGCATGGC
TGTTCTGTGCAACTGAGCGGCCGCACTGAGGTGGCTTATGTCTCCAAGGACACTCCTATGTTGCTTTACC
TCAACACTCACACAGCCTTGGAACAGATGCGGAGGCAAGCGGAAAAGGAAGAGCGAGGTCCCGGAGT
GATGGTAGTGGGCCCACTGATGTGGGCAAGTCTACAGTGTGTCGCCTTCTGCTCAACTACGCAGTGCCT
TTGGGCCGCGCTCCCACTTATGTGGAGCTGGATGTGGGCCAGGGTTCTGTGTCCATCCCTGGTACCATGG
GGGCCCTCTACATCGAGCGGCCTGCAGATGTGAAGAGGGTTTCTCTATCCAGGCCCTCTGGTGTATCA
TTTTGGTTCCCACTCCTGGCACTAACATCAAGCTTTATAAAGATTACATCTCGTTTAGCAGATGTG
TTCAACCAAAGGTGTGAGGTGAACCGAAGGGCATCTGTGAGTGGCTGTGTCATTAAACACCTGTGGCTGG
TCAAGGGCTCTGGTTACCAGGCTCTGGTGCATGCAGCCTCAGCTTTTGGGTGGATGTCGTTGTTGTTCT
GGATCAAGAACGACTGTACAATGAAGTAAACGGGACCTCCCCACTTTGTACGCACTGTGCTGCTCCCT
AAATCTGGGGGTGTGGTGGAGCGCTCCAAGGACTTCCGGCGGGAATGTAGGGATGAGCGTATCCGTGAGT
ATTTTTATGGATTCCGAGGCTGTTTCTATCCCATGCCTTCAATGTCAAATTTTCAAGATGTGAAATCTA
CAAAGTTGGGCAACCCACCATCCAGACTCCTGTTTACCTTTGGGCATGTCTCAAGAGGATAATCAGCTC
AAGCTAGTACCTGTCACTCCTGGGCGAGATATGGTGCACCACTACTGAGTGTAGCACTGCCGAGGGTA
CAGAGGAGAACCTGTCCGAGACAAGTGTAGCTGGCTTCATTGTGGTGACCACTGTGGACCTGGAGCATCA
GGTGTTTACTGTTCTGTCTCCAGCCCTCGCCCACTGCCTAAGAACTTCTTCTCATCATGGATATCCGG
TTCATGGATCTGAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201772 protein sequence Red=Cloning site Green=Tags(s) MGEEANDDKPPTTKFELERETELRFEVEASQSVQLELLTGMAEIFGTELTRNKKFTFDAGAKVAVFTWHG CSVQLSGRTEVAYVSKDTPMLLYLNHTALEQMRRQAEKEEERGPRVMVVGPTDVGKSTVCRLLLNYAVR LGRRPTYVELDVGGQSVSIPGTMGALYIERPADVEEGFSIQAPLVYHFGSTTPGTNIKLYNKITSRLADV FNQRCVENRRASVSGCVINTCGWVKSGSYQALVHAASAFEVDVVVLDQERLYNELKRDLPHFVRTVLLP KSGGVVERSKDFRRECRDERIREYFYGFRGCFYPHAFNVKFSVKIYKVGAPTIPDSCLPLGMSQEDNQL KLVPTPGRDMVHLLSVSTAEGTEENLSETSVAGFIVVTSVDLEHQVFTVLSPAPRPLPKNFLLIMDIR FMDLK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006831
ORF Size:	1275 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006831.1</u>
RefSeq Size:	1864 bp
RefSeq ORF:	1278 bp
Locus ID:	10978
UniProt ID:	<u>Q92989</u>
Cytogenetics:	11q12.1
MW:	47.6 kDa