

## Product datasheet for **RC201783L4V**

### MAP126 (SPAG5) (NM\_006461) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	MAP126 (SPAG5) (NM_006461) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MAP126
Synonyms:	DEEPEST; hMAP126; MAP126
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006461
ORF Size:	3579 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201783).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006461.3</a>
RefSeq Size:	3824 bp
RefSeq ORF:	3582 bp
Locus ID:	10615
UniProt ID:	<a href="#">Q96R06</a>
Cytogenetics:	17q11.2
Protein Families:	Druggable Genome, Stem cell - Pluripotency
MW:	134.4 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a protein associated with the mitotic spindle apparatus. The encoded protein may be involved in the functional and dynamic regulation of mitotic spindles. [provided by RefSeq, Jul 2008]