

## Product datasheet for **RC201785L1V**

### **FEN1 (NM\_004111) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	FEN1 (NM_004111) Human Tagged ORF Clone Lentiviral Particle
Symbol:	FEN1
Synonyms:	FEN-1; MF1; RAD2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_004111
ORF Size:	1140 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201785).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004111.4</a>
RefSeq Size:	2308 bp
RefSeq ORF:	1143 bp
Locus ID:	2237
UniProt ID:	<a href="#">P39748</a>
Cytogenetics:	11q12.2
Domains:	HhH2, XPG_N, XPG_I
Protein Families:	Druggable Genome, Stem cell - Pluripotency



[View online »](#)

**Protein Pathways:** Base excision repair, DNA replication, Non-homologous end-joining

**MW:** 42.6 kDa

**Gene Summary:** The protein encoded by this gene removes 5' overhanging flaps in DNA repair and processes the 5' ends of Okazaki fragments in lagging strand DNA synthesis. Direct physical interaction between this protein and AP endonuclease 1 during long-patch base excision repair provides coordinated loading of the proteins onto the substrate, thus passing the substrate from one enzyme to another. The protein is a member of the XPG/RAD2 endonuclease family and is one of ten proteins essential for cell-free DNA replication. DNA secondary structure can inhibit flap processing at certain trinucleotide repeats in a length-dependent manner by concealing the 5' end of the flap that is necessary for both binding and cleavage by the protein encoded by this gene. Therefore, secondary structure can deter the protective function of this protein, leading to site-specific trinucleotide expansions. [provided by RefSeq, Jul 2008]