

Product datasheet for **RC201787L3V**

DPAGT1 (NM_001382) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	DPAGT1 (NM_001382) Human Tagged ORF Clone Lentiviral Particle
Symbol:	DPAGT1
Synonyms:	ALG7; CDG-lj; CDG1j; CMS13; CMSTA2; D11S366; DGPT; DPAGT; DPAGT2; G1PT; GPT; UAGT; UGAT
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001382
ORF Size:	1224 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201787).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001382.2
RefSeq Size:	2150 bp
RefSeq ORF:	1227 bp
Locus ID:	1798
UniProt ID:	Q9H3H5
Cytogenetics:	11q23.3
Domains:	Glycos_transf_4



[View online »](#)

Protein Families:	Transmembrane
Protein Pathways:	Metabolic pathways, N-Glycan biosynthesis
MW:	46.1 kDa
Gene Summary:	The protein encoded by this gene is an enzyme that catalyzes the first step in the dolichol-linked oligosaccharide pathway for glycoprotein biosynthesis. This enzyme belongs to the glycosyltransferase family 4. This protein is an integral membrane protein of the endoplasmic reticulum. The congenital disorder of glycosylation type Ij is caused by mutation in the gene encoding this enzyme. [provided by RefSeq, Jul 2008]