

Product datasheet for RC201792A1V

Human UBE2N (NM_003348) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human UBE2N (NM_003348) AAV Particle
Tag: Myc-DDK
Symbol: UBE2N
Synonyms: HEL-S-71; UBC13; UbcH-ben; UbcH13; UBCHBEN
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC201792 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGC**

ATGGCCGGGCTGCCCCGAGGATCATCAAGGAAACCCAGCGTTTGCTGGCAGAACCAAGTTCCTGGCATCA
 AAGCCGAACAGATGAGAGCAACGCCGTTATTTTCATGTGGTCATTGCTGGCCCTCAGGATCCCCCTT
 TGAGGGAGGGACTTTTAACTTGAAGTATTCCTCCAGAAGAATACCAATGGCAGCCCTAAAGTACGT
 TTCATGACCAAAATTTATCATCCTAATGTAGACAAGTTGGGAAGAATATGTTTAGATATTTGAAAGATA
 AGTGGTCCCAGCACTGCAGATCCGCACAGTTCTGCTATCGATCCAGGCCTTGTTAAGTGCTCCCAATCC
 AGATGATCCATTAGCAAATGATGTAGCGGAGCAGTGGAAGACCAACGAAGCCCAAGCCATAGAAACAGCT
 AGAGCATGGACTAGGCTATATGCCATGAATAATATT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC201792 protein sequence
 Red=Cloning site Green=Tags(s)

MAGLPRIIKETQRLLAEPVPGIKAEPDESNAFYFHVVIAGPQDSPFEGGTFLKLEFLPEEYPMAPKVR
 FMTKIYHPNVDKLGRIKLDILKDKWSPALQIRTVLLSIQALLSAPNPDDPLANDVAEQWKTNEAQAIETA
 RAWTRLAMNNI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2


[View online »](#)

ACCN:	NM_003348
ORF Size:	456 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_003348.3</u>
RefSeq Size:	2568 bp
RefSeq ORF:	459 bp
Locus ID:	7334
UniProt ID:	<u>P61088</u>
Cytogenetics:	12q22
MW:	17.1 kDa