

## Product datasheet for **RC201792L1V**

### UBE2N (NM\_003348) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	UBE2N (NM_003348) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UBE2N
Synonyms:	HEL-S-71; UBC13; Ubch-ben; Ubch13; UBCHBEN
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_003348
ORF Size:	456 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201792).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003348.3</a>
RefSeq Size:	2568 bp
RefSeq ORF:	459 bp
Locus ID:	7334
UniProt ID:	<a href="#">P61088</a>
Cytogenetics:	12q22
Domains:	UBCc
Protein Families:	Druggable Genome



[View online »](#)

**Protein Pathways:** Ubiquitin mediated proteolysis

**MW:** 17.1 kDa

**Gene Summary:** The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. Studies in mouse suggest that this protein plays a role in DNA postreplication repair. [provided by RefSeq, Jul 2008]