

## Product datasheet for **RC201795L4V**

### PLK1 (NM\_005030) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PLK1 (NM_005030) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PLK1
Synonyms:	PLK; STPK13
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_005030
ORF Size:	1809 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC201795).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005030.3</a>
RefSeq Size:	2204 bp
RefSeq ORF:	1812 bp
Locus ID:	5347
UniProt ID:	<a href="#">P53350</a>
Cytogenetics:	16p12.2
Domains:	pkinese, POLO_box, TyrKc, S_TKc
Protein Families:	Druggable Genome, Protein Kinase




[View online »](#)

**Protein Pathways:** Cell cycle, Oocyte meiosis, Progesterone-mediated oocyte maturation

**MW:** 68.1 kDa

**Gene Summary:** The Ser/Thr protein kinase encoded by this gene belongs to the CDC5/Polo subfamily. It is highly expressed during mitosis and elevated levels are found in many different types of cancer. Depletion of this protein in cancer cells dramatically inhibited cell proliferation and induced apoptosis; hence, it is a target for cancer therapy. [provided by RefSeq, Sep 2015]

### Product images:



[RC201795L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC201795L4V particle to overexpress human PLK1-mGFP fusion protein.