

Product datasheet for **RC201796A1V**

Human HNRNPD (NM_031370) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human HNRNPD (NM_031370) AAV Particle
Tag:	Myc-DDK
Symbol:	HNRNPD
Synonyms:	AUF1; AUF1A; hnRNP0; HNRPD; P37
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201796 representing NM_031370 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGTCGGAGGAGCAGTTCGGCGGGACGGGCGGCGGCAGCGGCAACGGCGGCGGTAGGCGGCTCGGCGG
GCGAGCAGGAGGGAGCCATGGTGGCGGCACACAGGGGCGAGCGGCGGCGGGAAGCGGAGCCGGGAC
CGGGGGCGGAACCGCTCTGGAGGCACCGAAGGGGCGAGCGCCGAGTCGGAGGGGGCGAAGATTGACGCC
AGTAAGAACGAGGAGGATGAAGGCCATTCAAACCTCTCCCCACGACACTCTGAAGCAGCGACGGCACAGC
GGGAAGAATGAAAAATGTTTATAGGAGCCTTAGCTGGGACACTACAAAGAAAGATCTGAAGGACTACTT
TTCCAAATTTGGTGAAGTTGTAGACTGCACTCTGAAGTTAGATCCTATCACAGGGCGATCAAGGGGTTTT
GGCTTTGTGCTATTTAAGAATCGGAGAGTGTAGATAAGGTCATGGATCAAAAAGAACATAAATTGAATG
GGAAGGTGATTGATCCTAAAAGGGCCAAAGCCATGAAAACAAAAGAGCCGTTAAAAAATTTTTGTTGG
TGGCCTTTCTCCAGATACACCTGAAGAGAAAAAAGGGAGTACTTTGGTGGTTTTGGTGAGGTGGAATCC
ATAGAGCTCCCATGGACAACAAGACCAATAAGAGGCGTGGGTCTGCTTTATTACCTTTAAGGAAGAAG
AACCAGTGAAGAAGATAATGAAAAGAAATACCACAATGTTGGTCTTAGTAAATGTGAAATAAAGTAGC
CATGTCAAGGAACAATATCAGCAACAGCAACAGTGGGATCTAGAGGAGGATTTGCAGGAAGAGCTCGT
GGAAGAGGTGGTGGCCCCAGTCAAACCTGGAACAGGGATATAGTAACCTATTGGAATCAAGGCTATGGCA
ACTATGGATATAACAGCCAAGGTTACGGTGGTTATGGAGGATATGACTACACTGGTTACAACAACACTA
TGGATATGGTGATTATAGCAACCAGCAGAGTGTTATGGGAAGGTATCCAGGCGAGGTGGTCATCAAAAT
AGCTACAAACCATAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201796 representing NM_031370 Red=Cloning site Green=Tags(s) MSEEQFGGDGAAAAATAAVGGSAGEQEGAMVAATQGAAAAAGSGAGTGGGTASGGTEGGSASEGAKIDA SKNEEDEGHSNSSPRHSEAATAQREEWKMFIGGLSWDTTKKDLKDYFSKFGEVVDCTLKLDPIITGRSRGF GFVLFKESESVDKVMQKEHKLNGKVIDPKRAKAMKTKEPVKKIFVGGLSPDTPEEKIREYFGGFGEVES IELPMDNKTNKRGRFCFITFKEEEPVKKIMEKKYHNVGLSKCEIKVAMSKEQYQQQQWGSRGGFAGRAR GRGGGPSQNWQGYSNYWNQGYGNYGYNQGYGGYGGYDYGYNYYGYDYSNQSGYGKVSRRGGHQNSYKPY TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_031370
ORF Size:	1065 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_031370.2
RefSeq Size:	2257 bp
RefSeq ORF:	1068 bp
Locus ID:	3184
UniProt ID:	Q14103
Cytogenetics:	4q21.22
MW:	38.3 kDa