

Product datasheet for **RC201798A1V**

Human Proteasome beta 1 (PSMB1) (NM_002793) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Proteasome beta 1 (PSMB1) (NM_002793) AAV Particle
Tag:	Myc-DDK
Symbol:	Proteasome beta 1
Synonyms:	HC5; PMSB1; PSC5
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC201798 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTTGTCTCTACAGCCATGTATTCGGCTGCTGGCAGAGACTTGGGGATGGAACCGCACAGAGCCGCGG
GCCCTTTCAGCTGCGATTTTCGCCCTACGTTTTCAACGGAGGTACTATACTGGCAATTGCTGGAGAAGA
TTTTGCAATTGTTGCTTCTGATACTCGATTGAGTGAAGGGTTTTCAATTCATACGCGGATAGCCCCAAA
TGTTACAAATTAACAGACAAAACAGTCATTGGATGCAGCGTTTTTCATGGAGACTGTCTTACGCTGACAA
AGATTATTGAAGCAAGACTAAAGATGTATAAGCATTCCAATAATAAGGCCATGACTACGGGGCAATTGC
TGCAATGCTGTCTACAATCCTGTATTCAAGGCGCTTCTTCCATACTATGTTTACAACATCATCGGTGGA
CTTGATGAAGAAGGAAAGGGGGCTGTATACAGCTTTGATCCAGTAGGGTCTTACCAGAGAGACTCCTTCA
AGGCTGGAGGCTCAGCAAGTGCCATGCTACAGCCCCTGCTTGACAACCAGGTTGGTTTTAAGAACATGCA
GAATGTGGAGCATGTTCCGCTGTCCTTGGACAGAGCCATGCGGCTGGTGAAGATGTCTTCATTCTGCG
GCTGAGAGAGATGTGTACACTGGGACGCACTCCGGATCTGCATAGTGACCAAGAGGGCATCAGGGAGG
AACTGTTTCCTTAAGGAAGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC201798 protein sequence Red=Cloning site Green=Tags(s) MLSSTAMYSAGRDLMGEPHRAAGPLQLRFSPYVFNGGTILAIAGEDFAIVASDTRLSEGFSIHTRDSPK CYKLTDKTVIGCSGFHGDCLTLTKIIEARLKMYSNNKAMTTGAIAAMLSTILYSRRFFPYVYNIIGG LDEEGKGAIVSYFDPVGSYQRDSFKAGGSASAMLQPLLDNQVGFKNMQNVEHVPLSLDRAMRLVKDVFISA AERDVYTGDAALRICIVTKEGIREETVSLRKD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_002793
ORF Size:	723 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_002793.2</u>
RefSeq Size:	938 bp
RefSeq ORF:	726 bp
Locus ID:	5689
UniProt ID:	<u>P20618</u>
Cytogenetics:	6q27
MW:	26.5 kDa